

Digital Surgical Dentistry System

Revolutionizing Alveolar Height Augmentation:
Osteopore's Same-Day Customized Bone Regeneration Approach



Why Customize Your Osteomesh® for Alveolar Height Augmentation?

Traditionally, customizing a titanium mesh for alveolar height augmentation involves a lengthy process of 3D metal printing, requiring extended communication between the surgeon and the designer, and resulting in a delivery time of 6 – 9 weeks.

With Osteopore's Digital Surgical Dentistry (DSD) System, you now have the flexibility to **customize your Osteomesh® chair-side during surgery with precise shaping and moulding**, all within the same day as the patient's visit.

Streamline Surgical Procedures:

1 Same day treatment for patient, pre-moulding the Osteomesh® significantly **reduce the operation time**, benefiting both surgeons and patients.

Prioritize Soft Tissue Health:

2 Precise customization minimizes soft tissue disruption, ensuring **extra care for the soft tissue** by avoiding several trials of inserting the mesh.

Elevate Patient Comfort:

3 The Osteopore DSD System significantly reduces the duration of mouth opening, making the procedure more tolerable while **adding comfort to your patients**.

Are you ready to take your surgical dentistry practice to the next level?

Our Digital Surgical Dentistry (DSD) System combines a comprehensive package, including our state-of-the-art **Osteomesh®**, water bath, and the precision 3D3 Mini 2 printer. We offer a complete solution designed to enhance your surgical procedures and patient outcomes.

1 3D3 MINI 2 PRINTER:

Creating precise patient models at your fingertips

Engineered to perfection, this advanced printer is tailored for Alveolar Height Augmentation.

Equipped with impressive specifications, including:

- **4K Resolution:** Achieve anatomical accuracy with 4K clarity.
- **50µm per Layer:** Fine layers capture intricate features precisely.
- **Stability:** Double guide rail system ensures steady performance.
- **3 print platform sizes:**
 - i) 140mm*80mm
 - ii) 80mm*80mm
 - iii) 50mm*50mm
- **Print Height:** Accommodates complex anatomies with a 90mm height.
- **Light Cure Box:** Built-in for streamlined post-print curing.
- **Efficient Speed:** Optimal balance between speed and accuracy.
- **File Support:** Works seamlessly with STL files.
- **Multiple Connections:** Connect via LAN, WiFi, or USB.
- **Compact:** Dimensions of 261mm(W) x 322mm(D) x 420mm(H)
- **Weight:** 13kg
- **Warranty:** 1-year warranty





Enhancing Precision with Post Curing Station

Our DSD System includes a post curing station that complements the 3D3 Mini 2 printer.

Specifications:

- **Timer:** Tailor the curing time with intervals of 1 minute, ranging from 5 to 60 minutes.
- **Curing Method:** Choose between water and air for your specific needs.
- **Curing Wavelength:** 405nm
- **Curing Size:** 180X100mm
- **Dimensions:** 220mm (L) X 220mm (W) X 325mm (H)
- **Weight:** 5.5kg


2

PREPARATION OF OSTEOMESH®:

Revolutionizing Alveolar Height Augmentation

At the heart of our system lies the **Osteomesh®**— a semi-flexible, bioresorbable mesh engineered to provide exceptional dimensional retention and graft material support.

This innovative solution offers predictable shape and volume of regenerated bone tissue in guided bone regeneration, setting a new standard for surgical dentistry.



Porosity of Osteomesh®

3

MOULDING OF OSTEOMESH® ON WATER BATH BEFORE USE



Our water bath ensures the **Osteomesh®** achieves the ideal consistency, guaranteeing optimal results in Alveolar Height Augmentation.

It is recommended to place a sterile beaker or a sterile polyethylene bag/sheet in the water bath machine before pouring the sterile saline solution.

Customize Your Dental Scaffold for Alveolar Height Augmentation with Osteopore's Digital Surgical Dentistry (DSD) System now!

1 CRAFT PATIENT'S MODEL WITH 3D3 MINI 2 PRINTER



2 STERILIZATION OF THE MODEL

The model is printed from autoclavable polymer. After autoclaving the model, it can be used in the sterile field.



4 TRY-IN OF THE OSTEOMESH®

onto the patient's model crafted by the 3D3 Mini 2 Printer allowing for precise customization before the surgery.



3 PREPARATION OF OSTEOMESH®

- **SIZING OF DEFECT:** Select the suitable size of Osteomesh®, remove it from single wrap sterile packaging and place in dry sterile bowl.
- **TRIMMING OF OSTEOMESH®:** Using a pair of sterile surgical scissors, Osteomesh® is trimmed to the desired shape. Totally removing the borders of the mesh is recommended.
- **MOULDING OF OSTEOMESH®:** Warm up sterile saline to the appropriate temperature in the waterbath for moulding. Refer to the table for the appropriate temperature.

Mesh Thickness (mm)	Temperature (°C)
< 0.75	42 – 45
0.75 - 1.25	45 – 47
> 1.25	50

**Note: The saline temperature should not exceed the maximum temperature stated in the table.*

- ▶ Immerse Osteomesh® for approximately **10 seconds** to make it more malleable.
- ▶ Using the sterile model as a guide, contour and form the Osteomesh® in the desired vertical and horizontal measurements for approximately **5 seconds**.
- ▶ Remove Osteomesh® from the warm saline while holding it in its new shape for another **10 seconds**.
- ▶ Repeat step until desired shape is formed.




Scan the QR code for a video guide on trimming of the Osteomesh®.




5 INSERTING THE CUSTOMIZED, PRECISE OSTEOMESH® INTO THE PATIENT



PRODUCT SPECIFICATIONS For Ordering

Osteomesh®



	PRODUCT CODE	SIZE (L X B X T)/MM
	PC12(20,10,0.4)	20 x 10 x 0.4
	PC12(20,10,0.75)	20 x 10 x 0.75
	PC12(20,15,0.4)	20 x 15 x 0.4
	PC12(20,15,0.75)	20 x 15 x 0.75
	PC12(20,30,0.4)	20 x 30 x 0.4
	PC12(30,20,0.75)	30 x 20 x 0.75

Pictures of the product in actual size

Dental Model Material

PRODUCT CODE: SF001103B 1 Liter

- High Strength
- Autoclavable
- High Stability
- Grey, White, Tooth Color
- Water-washable



For professional use.

CAUTION: See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Osteopore®, the stylized Osteopore® logo, Osteoplug® are registered trademarks of Osteopore International Pte Ltd.

© 2023 Osteopore®. All rights reserved.

MCN 04 015B Rev 0

Osteopore®

Empowering Natural Tissue Regeneration

Osteopore International | SINGAPORE

2 Tukang Innovation Grove #09-06 JTC MedTech Hub, Singapore 618305